CE 1014

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CE 1014

Dental radiographic films DENTIX – instructions for use

Characteristic

DENTIX D is a double-emulsion, non-screen X-ray film of medium speed (ISO Class D), high contrast and fine grain, providing high image quality and outstanding detailed definition. The general purpose of the film is its use for the intraoral radiography in dental diagnostics.

DENTIX E is a double-emulsion high speed non-screen X-ray film (ISO Class E) with a high contrast and fine grain that provides high image quality with outstanding detail definition. The general purpose of the film is its use for the intraoral radiography in dental diagnostics.

DENTIX X-Stream is a double-emulsion extra high speed non-screen X-ray film (ISO Class F) with high contrast and low fine grain. The film has provided image quality and outstanding detail definition. The general purpose of the film is its use for the intraoral radiography in dental diagnostics.

Sides identification

Each film packet is protected with a lead sheet on the side opposite to the radiation source (marked on the back of the packet). Each film is provided with a protruding dot located near the edge of the film that serves as an identifier of the radiation side on the processed film. Its raised portion indicates the side facing the radiation source.

Darkroom illumination

The films are processed at yellow or diode safety illumination with a wavelength of 590 nm and higher. The length of exposure and a distance of the processed material from the illumination source should be tested.

Processing

DENTIX films can be processed manually or in processors. The FOMADENT concentrated solutions are recommended for processing the film to obtain the best results, however, any other processing chemicals of well-known brands can be used too. It is recommended to follow the producer's instructions while processing in processors. It is necessary to test the ideal setting of a specific device (the bath temperature, the movement speed etc.) in practice.

automatic processing	manual processing – developing
processing cycle/temperature	developing time/temperature
4,5–5,0 min / 25 °C	5,0 min / 20°C 4,0 min / 22°C 3,5 min / 24°C 1,5 min / 28°C
	Time of fixing/Temperature Not less than 2 min / 15–30°C

Note: It is necessary to wash the film in clean water at the end of manual processing.

Exposure conditions for DENTIX films (X-ray apparatus adjustment: 50–70 kV and 7–15 mA – use correct values recommended by the apparatus manufacturer)

Adjustment: 65 kV, 8 mA 20 cm focus-film distance				
	maxillary	exposure	mandibular	exposure
DENTIX D	incisor premolary molar	0,30 s 0,40 s 0,45 s	incisor premolary molar	0,25 s 0,28 s 0,30 s
DENTIX E	incisor premolary molar	0,18 s 0,23 s 0,25 s	incisor premolary molar	0,14 s 0,14 s 0,18 s
DENTIX X-Stream	incisor premolary molar	0,14 s 0,18 s 0,20 s	incisor premolary molar	0,11 s 0,12 s 0,14 s

For making exposures of children reduce the exposure time approx. by 33%.
For making exposures of empty patches reduce the exposure time approx. by 25%.
For obtaining the best possible results all the necessary changes of exposure parameters (i.e. exposure time, mA, kV or any changes of the focus-film distance) shall be reflected in other parameters.
On X-ray devices with a DC generator, to reduce exposure by about 33%.

Storage

DENTIX films should be stored in its original packaging in a cool dry place at a temperature from 10°C to 25°C and a relative humidity of 40–60 % protected from damaging fumes, gases and ionizing radiation. It is recommended to keep the film in a refrigerator with the temperatures between 5 and 8°C in case of a long-term storage. Before use, the film in the intact original packaging should be allowed to adjust to room conditions for at least 2–4 hours.

Warning

- The product is not sterile ; in accordance with the hygienic rules the disinfection of individual film packagings (each packet) must be applied before and after use and the procedure intended for specific healthcare department must be obeyed.
- Usage and processing of the DENTIX films result in waste that is classified as hazardous waste and for this reason an ecological liquidation and recycling is necessary.
Classification of the waste:
polyvinylchloride waste, waste of developers, lead waste, waste of fixing bath
- The manipulation with the waste should be in accordance with valid statutory instruments. It is necessary to work with individual packets carefully, not to flex them nor to push them down using nails or sharp tools. The film inside the packet is sensitive to pressure and the excessive pushing can cause pressure records which may complicate the picture evaluation.
- Do not use if the packaging is damaged.

Please apply any inquiries or reminders to
FOMA BOHEMIA spol. s r.o., Jana Krušinky 1737/6, 500 02 Hradec Králové, Czech Republic
phone: +420 495 733 320, e-mail: foma@foma.cz, www.foma.cz

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